

Identifying Causes of Steam Sterilization Process Failures

Key variables to consider when investigating a steam sterilization process failure

Human Errors

Incorrect use and interpretation of monitoring tools

- ☐ Not determining that the physical monitors were correct for the load
- ☐ Not identifying that the correct cycle was run for the load contents
- ☐ Using the wrong BI/PCD for the load
- ☐ Using the wrong Class 5 integrating indicator/PCD for the load
- ☐ Not following the BI/PCD or Class 5 integrating indicator instructions for use
- ☐ Incorrect reading of Class 5 integrating indicator or BI result
- ☐ Using the wrong internal chemical indicator for the cycle
- ☐ Incorrect reading of internal chemical indicator
- ☐ Incorrect storage of any chemical or biological indicators

Improper cycle for the load contents

- ☐ Not following the container or instrument/container manufacturers' instructions for use
- ☐ Not verifying instrument/container manufacturers' sterilization parameters for use in your sterilizers using CSA product testing protocol.

Inappropriate packaging materials or packaging technique

- ☐ Incorrect packing or container system for the cycle parameters
- ☐ Not correctly preparing the container for use (i.e., filters and valves or appropriate bottom tray)
- ☐ Using an oversized wrap for application
- ☐ Not using a mesh bottom perforated tray that allows air removal and steam penetration
- ☐ Placing a paper-plastic pouch inside a container system or wrapped tray

- ☐ Placing a folded paper-plastic pouch inside another paper-plastic pouch
- ☐ Preparing textile packs that are too dense to sterilize in the cycle parameters chosen
- ☐ Not placing basins in same direction
- ☐ Not using non-linting absorbent material between nested basins
- ☐ Not disassembling or opening hinged instruments or surgical supplies
- ☐ Not holding packaging materials at 68°F-73°F (20°C-23°C), 30-60% RH for 2 hours prior to use

Sterilizer loading

- ☐ Stacking container systems if not recommended by manufacturer
- ☐ Stacking perforated instrument trays
- ☐ Not laying instrument trays flat or parallel to the shelf
- ☐ Laying peel pouches flat instead of on edge; not properly spaced or with plastic sides not facing one direction
- ☐ Not placing basins on edge
- ☐ Not placing fabric packs on edge
- ☐ Placing packages too close to each other impeding air removal and sterilant penetration around and through load

Poor Steam Quality or Quantity

Wet steam

- ☐ Improperly insulated steam lines
- ☐ Malfunctioning or missing trap in steam line
- ☐ Malfunctioning or no drain check valve
- ☐ Steam contact with a cold load
- ☐ Too much water in steam produced at boiler

Continued on reverse

Poor Steam Quality or Quantity (cont.)

Superheated steam

- ☐ Improper chamber heat up
- ☐ Desiccated packaging materials (e.g., towels)
- ☐ Steam pressure too low for the temperature
- ☐ Excessive reduction of steam pressure too close to sterilizer
- ☐ Faulty steam control valve
- ☐ Faulty pressure reducer control valve

Other reasons

- ☐ Variations in steam pressure due to clogged filter, poorly engineered piping or excessive demands
- ☐ Pressure gauges and controllers out of calibration
- ☐ Clogged steam lines
- ☐ Clogged steam supply strainer
- ☐ Clogged chamber drain line, strainer or chamber drain screen
- ☐ Malfunctioning valves

Incomplete Air Removal

- ☐ Inadequate vacuum or vacuum depth or other air removal system
- ☐ Clogged chamber drain line, strainer or chamber drain screen
- ☐ Clogged vent lines
- ☐ Leak caused by a faulty door gasket
- ☐ Leak in other areas of chamber
- ☐ Plugged, faulty or maladjusted control valves

- ☐ Low steam pressure
- ☐ Too high water temperature
- ☐ Inadequate water supply pressure
- ☐ Clogged water supply strainer
- ☐ Air trapped by the load
- ☐ Incorrect cycle parameters for the load

Inadequate Cycle Temperature

- ☐ Temperature gauge out of calibration
- ☐ Long heat-up time of large loads (i.e., heat lag)
- ☐ Clogged chamber drain line, strainer or chamber drain screen
- ☐ Variations in steam pressure due to clogged filter, poorly engineered piping or excessive demands on the steam supply
- ☐ Presence of non-condensable gases in steam line and load
- ☐ Inadequate steam supply pressure
- ☐ Clogged steam supply strainer

Insufficient Time at Temperature

- ☐ Control timer out of calibration
- ☐ Inappropriate cycle parameters for the load being processed
- ☐ Come up time less than 1.5 minutes in a gravity 270°F to 275°F (132°C to 135°C) 3 minute cycle
- ☐ Oversized load

**3M STERILIZATION
ASSURANCE**

**EVERY STEP. EVERY DETAIL.
EVERY DAY.**



3M Health Care
3M Canada Company
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1 800-364-3577
www.3m.com/ca/healthcare

3M Infection Prevention Division
3M Health Care
3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
U.S.A.
1 800 228-3957
www.3m.com/healthcare

Attest and 3M are trademarks of 3M.
Used under license in Canada.
Please recycle. Printed in Canada
© 3M 2009. All rights reserved.
0904-5502E